

Progesterone Receptor antibody [1A6]

Cat. No. GTX75363

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P
Reactivity	Human

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/25-1/50
IHC-P	1/40

Note : This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Not tested in other applications.

Calculated MW 99 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide corresponding to a site of predicted high antigenicity on the human progesterone receptor.
Purification	Tissue Culture Supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX75363 IHC-P Image


IHC-P analysis of breast cancer biopsy tissue using GTX75363 Progesterone Receptor antibody. Heat mediated antigen retrieval using citrate buffer pH6.0.


GTX75363 IHC-P Image


IHC-P analysis of breast cancer biopsy tissue using GTX75363 Progesterone Receptor antibody. Heat mediated antigen retrieval using citrate buffer pH6.0.


GTX75363 IHC-P Image

IHC-P analysis of breast cancer biopsy using GTX75363 Progesterone Receptor antibody [1A6].


GTX75363 IHC-P Image

IHC-P analysis of breast cancer biopsy tissue using GTX75363 Progesterone Receptor antibody. Heat mediated antigen retrieval using citrate buffer pH6.0.


GTX75363 IHC-P Image

IHC-P analysis of breast cancer biopsy tissue using GTX75363 Progesterone Receptor antibody. Heat mediated antigen retrieval using citrate buffer pH6.0.



For full product information, images and publications, please visit our [website](https://www.genetex.com).