

CD81 antibody [Eat 2] (FITC)

Cat. No. GTX75431**Host** Armenian Hamster**Clonality** Monoclonal**Isotype** IgG**Applications** FCM**Reactivity** Mouse**Package**

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

FCM

Recommended dilution

Assay dependent

Not tested in other applications.


Properties

Form Liquid**Buffer** PBS, BSA**Preservative** 0.02% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.**Concentration** 0.1 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** 38C13, murine B cell line**Purification** Purified IgG**Conjugation** Fluorescein isothiocyanate (FITC) [Wavelength](#)**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX75431 FCM Image

Staining of mouse spleen cells with Hamster anti Mouse CD81:FITC (GTX75431)



For full product information, images and publications, please visit our [website](#).