

## Mouse Anti-Rat IgE (Heavy chain) antibody [MARE-1] (HRP)

## Cat. No. GTX75595

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, ELISA
Reactivity	Rat

References ( 2 ) Package 250 µg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	0.5μg/ml

Not tested in other applications.

Properties	
Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rat IR162, IR1016, IR2 and IR410 IgE myeloma proteins.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Horseradish peroxidase(HRP)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 03 Page 1 of 1