

# CD200 antibody [MRC OX90]

**Cat. No. GTX75633**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	IHC-Fr, FCM, ELISA
<b>Reactivity</b>	Mouse

**Package**  
250 µg

## Applications

### Application Note

For ELISA: Use at an assay dependent dilution. For FACS: Use at an assay dependent dilution. For IHC-Fr: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Mouse CD200-rat CD4 fusion protein
<b>Purification</b>	Protein G purified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).