

CD173 antibody [BRIC231]

Cat. No. GTX75683

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM, Haemagglutination
Reactivity	Human, Pig

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
Haemagglutination	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris-Glycine
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human erythroleukemic cell line (HEL) established from a 30 year old patient with relapsed erythroleukemia following treatment for Hodgkin lymphoma.
Purification	Purified IgG
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).