

p15 INK4b antibody [DCS114.1]

Cat No. GTX75687

| | |
|--------------------|---------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Application | WB, IHC-Fr, IP, IHC |
| Reactivity | Human |

Reference (2)
 Package
 50 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Dilution |
|--------------------|-----------------|
| WB | Assay dependent |
| IHC-Fr | Assay dependent |
| IP | Assay dependent |
| IHC | Assay dependent |

Not tested in other applications.

Calculated MW 15 kDa. ([Note](#))

PROPERTIES

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS, 0.08% Sodium Azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1.0 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Full length recombinant p15 His tagged |
| Purification | Purified IgG |
| Conjugation | Unconjugated |
| Note | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |



For full product information, images and publications, please visit our [website](#).