

p15 INK4b antibody [DCS114.1]

Cat. No. GTX75687

Host	Mouse	
Clonality	Monoclonal	
Isotype	lgG1	
Applications	WB, IHC-Fr, IP, IHC	
Reactivity	Human	

References (2) Package 50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 15 kDa. (Note)

Properties		
Form	Liquid	
Buffer	PBS	
Preservative	0.08% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)	
Immunogen	Full length recombinant p15 His tagged	
Purification	Purified IgG	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 02 Page 1 of 2



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 02 Page 2 of 2