

HLA-G antibody [MEM-G/1]

Cat. No. GTX75708

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC
Reactivity	Human

Package

200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Product Note Recognises denatured free HLA-G but does not recognise HLA-G bound to the cell surface.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Denatured heavy chain full-length recombinant human HLA-G
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX75708 Image

The antibody stains infiltrating macrophages in pulmonary diseases. In the top left corner see detail of macrophage.



For full product information, images and publications, please visit our [website](#).