

Integrin alpha 4 antibody [R1-2]

Cat. No. GTX75782

Host	Rat
Clonality	Monoclonal
Isotype	IgG2b
Applications	FCM, IP, IHC
Reactivity	Mouse

References (1)

Package

250 µg

Applications

Application Note

For FACS: Use at an assay dependent dilution. For IHC: Use at an assay dependent dilution. For IP: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peyers Patch HEV binding lymphoma line (TK1)
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).