

CD41 antibody [MWReg30]

Cat. No. GTX76011

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P, IHC-Fr, FCM, IP, IHC
Reactivity	Human, Mouse

References (18)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	1/50-1/100
IP	Assay dependent
IHC	Assay dependent

Note : Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified murine platelets.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX76011 FCM Image

FACS analysis of mouse peripheral blood platelets using GTX76011 CD41 antibody [MWReg30].



GTX76011 IHC-Fr Image

The data was published in the journal PLoS One in 2013. [PMID: 23991201](https://pubmed.ncbi.nlm.nih.gov/23991201/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).