

## Mouse Anti-Rat IgG2a (Heavy chain) antibody [MARG2a-1] (HRP)

Cat. No. GTX76070

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	ELISA
Reactivity	Rat

Package  
250 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	0.5µg/ml

Not tested in other applications.

## PROPERTIES

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rat IgG2a
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Horseradish peroxidase(HRP)

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).