

CD45RA antibody [OX-33] (FITC)

Cat. No. GTX76085

| | |
|--------------|------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | FCM |
| Reactivity | Rat |

Package
100 µg

Applications

Application Note

For FACS: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Product Note Rat CD45RA/B

Properties

| | |
|---------------|--|
| Form | Liquid |
| Buffer | PBS, BSA |
| Preservative | 0.02% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light. |
| Concentration | 0.1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Purified Rat spleen L-CA |
| Purification | Protein G purified |
| Conjugation | Fluorescein isothiocyanate (FITC) Wavelength |

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx76085 FCM Image

Staining of rat spleen cells with Mouse anti Rat CD45RA: FITC (GTx43573)



For full product information, images and publications, please visit our [website](#).