

CD45RA antibody [OX-33] (PE)

Cat. No. GTX76087

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	FCM
Reactivity	Rat

Package 50 μg

Applications

Application Note

For FACS: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Rat CD45RA/B **Product Note**


Properties	
Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified rat LCA
Purification	Purified IgG
Conjugation	Phycoerythrin (PE) <u>Wavelength</u>
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 28 Page 1 of 2

DATA IMAGES



GTX76087 FCM Image

Staining of rat spleen cells with Mouse anti Rat CD45RA: FITC (GTX43573)



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 10 / 28 Page 2 of 2