

CD45 antibody [OX-1] (Biotin)

Cat. No. GTX76132

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	IHC-Fr, FCM
Reactivity	Rat

References (2) Package 1 ml

Applications

Application Note

For FACS: Use at an assay dependent dilution. Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the researcher.

Rat Leukocyte Common Antigen **Product Note**

Properties	
Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Thymocyte membrane proteins
Purification	Protein G purified
Conjugation	Biotin
Conjugation	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 28 Page 1 of 1