

CD5 antibody [MRC OX-19] (PE)

Cat. No. GTX76270

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS
Reactivity	Rat

Package

50 µg

APPLICATION

Application Note

For FACS: Use 0.5 µg to label 10⁶ cells. Protocols are available upon request. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	53 kDa. (Note)
Product Note	Rat CD5

PROPERTIES

Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rat thymocyte glycoproteins
Purification	Protein G affinity purified
Conjugation	Phycoerythrin (PE)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).