

## CD45RC antibody [OX-22] (PE)

Cat. No. GTX76278

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM
Reactivity	Rat

## Package

50 µg

## Applications

## Application Note

For FACS: Use 0.1 µg to label 10<sup>6</sup> cells. Optimal dilutions/concentrations should be determined by the end user.

**Product Note** rat CD45RC

## Properties

Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Phytohaemagglutinin (PHA) -activated rat lymphocytes
Purification	Protein G affinity purified
Conjugation	Phycoerythrin (PE) <a href="#">Wavelength</a>

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).