

CD43 antibody [W3/13] (PE)

Cat. No. GTX76288

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS
Reactivity	Rat

Package
50 µg

APPLICATION

Application Note

For FACS: Use 0.5 µg to label 10⁶ cells. Protocols are available upon request. Optimal dilutions/concentrations should be determined by the end user.

Product Note Rat CD43

PROPERTIES

Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rat thymocyte membrane
Purification	Protein G affinity purified
Conjugation	Phycoerythrin (PE)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX76288 FACS Image

Staining of rat splenocytes with Mouse anti Rat CD43: FITC (GTX76286)



For full product information, images and publications, please visit our [website](#).