

## CD43 antibody [W3/13] (PE)

**Cat. No. GTX76288**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	FCM
<b>Reactivity</b>	Rat

**Package**  
50 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	0.5 µg for 10 <sup>6</sup> cells

Not tested in other applications.

<b>Product Note</b>	Rat CD43
---------------------	----------

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Rat thymocyte membrane
<b>Purification</b>	Protein G affinity purified
<b>Conjugation</b>	Phycoerythrin (PE) <a href="#">Wavelength</a>

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX76288 FCM Image**

Staining of rat splenocytes with Mouse anti Rat CD43: FITC (GTX76286)



For full product information, images and publications, please visit our [website](#).