

## CD4 antibody [W3/25] (APC)

Cat. No. GTX76315

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS
Reactivity	Rat

Reference ( 3 )

Package

50 test

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	1/5-1/10

**Note : Use 10 $\mu$ l of the suggested working dilution to label 10<sup>6</sup> cells in 100 $\mu$ l.**

Not tested in other applications.

**Calculated MW** 51 kDa. ( [Note](#) )

## PROPERTIES

Form	Liquid
Buffer	PBS, 1% BSA, 5% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Rat Thymocyte Membrane Glycoproteins
Purification	Purified IgG
Conjugation	Allophycocyanin (APC)

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).