

Mouse IgG2b isotype control (Biotin)

Cat. No. GTX76398

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	FCM

Package
100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

FCM

Recommended dilution

Assay dependent

Not tested in other applications.

Product Note

This antibody is suitable as a negative control with human, bovine, ovine, porcine, canine and guinea-pig tissues. And it is not suitable for use as a negative control on rat tissues.

Properties

Form Liquid**Buffer** PBS, 1% BSA**Preservative** 0.09% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.1 mg/ml (Please refer to the vial label for the specific concentration.)**Purification** Protein A purified
From tissue culture supernatant**Conjugation** Biotin**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).