

Mouse IgG2b isotype control (Azide free)

Cat. No. GTX76402

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	IHC-P, IHC-Fr, FCM, IHC (Resin sections), RIA

Package
250 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IHC (Resin sections)	Assay dependent
RIA	Assay dependent

Not tested in other applications.

Product Note

This antibody is suitable as a negative control with human, bovine, ovine, porcine, canine and guinea-pig tissues. And It is not suitable for use as a negative control on rat tissues.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).