

CD45RA antibody [F8-11-13] (PE)

Cat. No. GTX76598

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	FCM
Reactivity	Human, Rhesus Monkey

Package 100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher

Optimal anations/co	The chitations should be determined by the rescarcion.		
Suggested dilution	Recommended dilution		
FCM	Neat		
Note : Use 10μl of the suggested working dilution to label 10 ⁶ cells or 100μl whole blood.			
Not tested in other ap	uplications.		
Product Note	This antibody reacts with the high molecular weight form of the leucocyte common antigen (LCA) expressed by B		

Product Note	lymphocytes and a subset of T lymphocytes.

Properties	
Form	Liquid
Buffer	PBS, 1% BSA, 5% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human T lymphocytes
Purification	Protein G purified From tissue culture supernatant
Conjugation	Phycoerythrin (PE) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 28 Page 1 of 1