

Mouse IgG1 isotype control (PE)

Cat. No. GTX76639

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM

Package
50 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent

Not tested in other applications.

Product Note

This antibody is suitable as a negative control with human, bovine, ovine, porcine, equine, canine, lapine and guinea-pig tissues. And it is not suitable for use as a negative control on rat tissues.

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 5% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Purification	Protein A purified From tissue culture supernatant
Conjugation	Phycoerythrin (PE) Wavelength

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).