

Human IgG1 Lambda chain isotype control

Cat. No. GTX76676

Host	Human
Applications	WB, Dot, ELISA, Immunoassay, Immunodiffusion, Immunoelectrophoresis

References (1)
Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
Dot	Assay dependent
ELISA	Assay dependent
Immunoassay	Assay dependent
Immunodiffusion	Assay dependent
Immunoelectrophoresis	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	20mM TBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.15 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Fractionation/ion-exchange/affinity chromatography From human plasma
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).