

## Rabbit Anti-Bovine Whole Serum antibody

**Cat. No. GTX77157**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	Immunoelectrophoresis
<b>Reactivity</b>	Bovine

**Package**  
1 ml

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunoelectrophoresis	Assay dependent

Not tested in other applications.

**Product Note**

This antiserum forms multiple arcs against normal bovine serum by immunoelectrophoresis and has strong anti-albumin reactivity.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	Serum
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Immunogen</b>	Pooled normal bovine serum in Freund's adjuvant.
<b>Purification</b>	Unpurified
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).