

Rabbit Anti-Rat IgG antibody

Cat. No. GTX77340

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Reactivity	Rat

Package
500 µg

PROPERTIES

Form	Liquid
Buffer	50mM HEPES, 100mM NaCl, 0.2% BSA
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).