

Hemoglobin antibody

Cat. No. GTX77484

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

References (2)

Package

1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	2 ug/ml

Not tested in other applications.

Product Note GTX77484 forms a single precipitin band with human adult hemoglobin A1, A2, F and S, by immunoelectrophoresis.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Goats were immunized with human hemoglobin in Freund's adjuvant. After multiple immunizations the goats were bled, serum collected and the anti hemoglobin immunoaffinity purified.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX77484 ELISA Image

ELISA analysis of immunogen using GTX77484 Hemoglobin antibody.



For full product information, images and publications, please visit our [website](#).