

EGFR antibody

Cat. No. GTX77610

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2b
Applications	WB, ICC/IF, IHC-Fr, FCM, IP
Reactivity	Human, Horse

Package 100 μg

Applications

Application Note

Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 134 kDa. (Note)

Product Note This antibody does not react with erbB-2, erbB-3 and erbB-4.

Properties	
Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human epidermoid carcinoma line A431; epitope mapped between aa 6-273 of human EGFR
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com