

UBIAD1 antibody

Cat. No. GTX77653

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

Western blot - Use at 1:500 dilution. However, the investigator should determine the optimal dilution for a specific application such as immunohistochemistry.

Product Note Reacts with residues 32-45 [PEQDRLPQRSWRQK] of the 37 kDa human UBIAD1 protein. This sequence is 85% identical in mouse and rat.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide: PEQDRLPQRSWRQK
Purification	Affinity purified with antigen. From polyclonal serum.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).