

mtTFA antibody, N-term

Cat. No. GTX77852

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

References (2)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml

Not tested in other applications.

Calculated MW 29 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human mtTFA
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx77852 WB Image


WB analysis of Jurkat cells using GTx77852 mtTFA antibody at 0.2-1µg/ml.


GTx77852 WB Image

WB analysis of bEND.3 cells using GTx77852 mtTFA antibody at 1:1000.

Lane 1 : bEND.3 cells (10µg)

Lane 2 : bEND.3 cells (20µg)


GTx77852 WB Image

WB analysis of HepG2, Hep3B, and K562 cells using GTx77852 mtTFA antibody at 1:1000.

Lane :

1: HepG2 cell lysate

2: HepG2 cell lysate

3: Hep3B cell lysate

4: Hep3B cell lysate

5: K562 cell lysate

6: K562 cell lysate



For full product information, images and publications, please visit our [website](https://www.genetex.com).