

## MYF5 antibody, N-term

Cat. No. GTX77867

|              |              |
|--------------|--------------|
| Host         | Rabbit       |
| Clonality    | Polyclonal   |
| Isotype      | IgG          |
| Applications | WB, ICC/IF   |
| Reactivity   | Human, Mouse |

References ( 1 )

Package

50 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.2-2.5 ug/ml        |
| ICC/IF             | Assay dependent      |

Not tested in other applications.

Calculated MW 28 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 2% Sucrose                                                                                                                                                                                                            |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                               |
| Immunogen     | A synthetic peptide corresponding to a N-terminal region of Human MYF5                                                                                                                                                     |
| Purification  | Affinity Purified                                                                                                                                                                                                          |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX77867 WB Image

WB analysis of mouse skeletal and muscle tissue using GTX77867 MYF5 antibody at 1:1000.

Lane1: Mouse heart lysate

Lane2: Mouse skeletal muscle lysate



### GTX77867 WB Image

WB analysis of Jurkat cells using GTX77867 MYF5 antibody at 1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).