

FOXE3 antibody, C-term

Cat. No. GTX77875

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	2-10 µg/ml

Not tested in other applications.

Product Note We do not recommend use of this product for Zebrafish samples.

Properties


Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a C-terminal region of Human FOXE3
Purification	Affinity Purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX77875 IHC-P Image

IHC-P analysis of human muscle tissue using GTX77875 FOXE3 antibody at 4.0-8.0μg/ml.


GTX77875 IHC-P Image

IHC-P analysis of human spermatophore tissue using GTX77875 FOXE3 antibody at 4.0-8.0μg/ml.



For full product information, images and publications, please visit our [website](#).