

GPR23 antibody

Cat. No. GTX77987

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, FCM
Reactivity	Human, Mouse, Rat, Hamster, Monkey, Bat

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	2 µg/ml
FCM	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 19 amino acid peptide from N-terminal extracellular domain of human LPAR4 / GPR23.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX77987 IHC-P Image


IHC-P analysis of ovary, carcinoma tissue using GTX77987 GPR23 antibody.


GTX77987 IHC-P Image

IHC-P analysis of human ovary, carcinoma tissue using GTX77987 GPR23 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX77987 IHC-P Image

IHC-P analysis of human brain neuron tissue using GTX77987 GPR23 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX77987 IHC-P Image

IHC-P analysis of human kidney tissue using GTX77987 GPR23 antibody.
Antigen retrieval : Heat-induced antigen retrieval



For full product information, images and publications, please visit our [website](#).