

5-HT7 receptor antibody

Cat. No. GTX78121

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	13 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 16 amino acid peptide from C-terminus of human 5HT7 Receptor.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES


GTx78121 IHC-P Image

IHC-P analysis of human brain, cortex, neurons tissue using GTx78121 5-HT7 receptor antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 19 µg/ml


GTx78121 IHC-P Image

IHC-P analysis of human skeletal muscle tissue using GTx78121 5-HT7 receptor antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 19 µg/ml



For full product information, images and publications, please visit our [website](#).