

CDC7 antibody

Cat. No. GTX78149

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (1)

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	5.2 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 19 amino acid peptide from kinase domain of human CDC7.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX78149 IHC-P Image


IHC-P analysis of human prostate, carcinoma tissue using GTX78149 CDC7 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX78149 IHC-P Image


IHC-P analysis of brain, alzheimer's neurofibrillary tangle and granulovacuolar degeneration tissue using GTX78149 CDC7 antibody.


GTX78149 IHC-P Image

IHC-P analysis of colon, adenocarcinoma tissue using GTX78149 CDC7 antibody.


GTX78149 IHC-P Image

IHC-P analysis of human ovary, carcinoma tissue using GTX78149 CDC7 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX78149 IHC-P Image

IHC-P analysis of human testis tissue using GTX78149 CDC7 antibody.
Antigen retrieval : Heat-induced antigen retrieval



For full product information, images and publications, please visit our [website](https://www.genetex.com).