

DPP10 antibody

Cat. No. GTX78173

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

Recommended dilution for Western blot 1:500-1:1,000. However, the investigator should determine the optimal dilution for a specific application such as immunohistochemistry.

Calculated MW	91 kDa. (Note)
Product Note	Reacts with residues 521-539 LKEAILKKKIGKPEIKILH of human DPP10. Other species not tested.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Affinity purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx78173 WB Image

Western Blot detection of Flag-DPP10 in transfected CHO cells with anti-DPP10 antibody at a 1:1000 dilution. The expected MW is 75 kDa; however, oligomerization and differential glycosylation produces other bands as observed with other DPP family members.



For full product information, images and publications, please visit our [website](https://www.genetex.com).