

JNK (phospho Thr183/Tyr185) antibody

Cat. No. GTX78197

Host	Chicken
Clonality	Polyclonal
Isotype	IgY
Applications	ICC/IF, ELISA, IHC
Reactivity	Human

Package
100 µl

Applications

Application Note

The product can be used at 1:1,000 to 1:4,000 for immunoblot testing and ELISA or 1:500 to 1:2,000 for UHC or ICC.

Product Note

This antibody is raised against JNK1/2 phosphorylated at Thr183/Tyr185. Based on sequence homology, it is predicted to react with JNK3 when phosphorylated at the corresponding residues - Thr221/Tyr223.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Phospho-JNK (pT183-P-pY185) peptide
Purification	Affinity purified with antigen.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).