

Rb2 / p130 antibody [DCS211.6]

Cat. No. GTX78225

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	WB
Reactivity	Human

Package
100 µg

APPLICATION

Application Note

Western blot (1.0 ug/ml)

Calculated MW 128 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	100 ug purified IgG1 with 50% glycerol, 0.01% sodium azide and 1.0 mg/ml BSA.
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Balb/c mice were immunized with human Rb2 / p130 protein and fusing the splenocytes with Sp2/0 mouse myeloma cells.
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).