

Rb2 / p130 antibody [DCS211.6]

Cat. No. GTX78225

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

Western blot (1.0 ug/ml)

Calculated MW 128 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	100 ug purified IgG1 with 50% glycerol, 0.01% sodium azide and 1.0 mg/ml BSA.
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Balb/c mice were immunized with human Rb2 / p130 protein and fusing the splenocytes with Sp2/0 mouse myeloma cells.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).