

Musashi 1 antibody

Cat. No. GTX78273

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC
Reactivity	Human, Mouse, Rat

References (2)
Package
100 µg

Applications

Application Note

Immunocytochemistry and Immunohistochemistry - Use at 1:100 dilution. Western blot - Use at 1:500 dilution. However, the investigator should determine the optimal dilution for a specific application such as immunohistochemistry.

Product Note Human, mouse and rat Musashi-1 protein

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Reacts with residues 5-21 [APQPGLASPDSPHPDCK] of the human, mouse and rat Musashi 1 protein.
Purification	Affinity purified IgG
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES


Image of staining with Musashi polyclonal antibody (1:100 dilution) shows Musashi (green) staining on neural rosettes (human). Nuclei are counterstained blue (DAPI).

GTX78273 ICC/IF Image


For full product information, images and publications, please visit our [website](#).