

CD46 antibody [MEM-258] (PE)

Cat. No. GTX78353

| | |
|---------------------|---------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | FCM |
| Reactivity | Human, Bovine |

References (2)

Package

100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|--|
| FCM | 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension |

Note : Do not need a permeabilization step.

Not tested in other applications.

Product Note This antibody recognizes an epitope on SCR4 (the membrane-proximal SCR) domain of CD46 (Membrane cofactor protein).

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 15mM Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | HPB-ALL human T cell line |
| Purification | Purified by size-exclusion chromatography |
| Conjugation | Phycoerythrin (PE) Wavelength |


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX78353 FCM Image**

FACS analysis of human peripheral blood using GTX78353 CD46 antibody [MEM-258] (PE).



For full product information, images and publications, please visit our [website](#).