

HLA-DR1 antibody [MEM-267]

Cat. No. GTX78393

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))

Product Note The antibody MEM-267 specifically binds to the empty but not peptide-loaded form of HLA-DR1.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified, insoluble DR1 beta chain (DRB1*0101) expressed in E. coli inclusion bodies.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX78393 WB Image

WB analysis of Raji (A) and Jurkat (B) cell lysates using GTX78393 HLA-DR1 antibody [MEM-267].



For full product information, images and publications, please visit our [website](https://www.genetex.com).