

HLA-G antibody [MEM-G/2]

Cat. No. GTX78394

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IHC-Fr
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent

Not tested in other applications.

Calculated MW	38 kDa. (Note)
Product Note	This antibody recognizes the free heavy chain of all the HLA-G isoforms.

Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Bacterially expressed recombinant HLA-G heavy chain (denatured).
Purification	Purified by precipitation and chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX78394 WB Image

WB analysis of HLA-G1 transfectants (LCL-HLA-G1) using GTX78394 HLA-G antibody [MEM-G/2].



GTX78394 IHC-P Image

IHC-P analysis of kidney allograft biopsy tissue using GTX78394 HLA-G antibody [MEM-G/2].



For full product information, images and publications, please visit our [website](#).