

Mouse IgM isotype control [PFR-03]

Cat. No. GTX78439

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	FCM

Package
100 µg

Applications

Application Note

Optimal dilutions/concentrations should be determined by the end user.

Product Note PFR-03 reacts with undefined epitope on a plant pathogen.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein G affinity purified. From hybridoma culture supernatant.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).