

Presenilin 1 antibody [34-4C3]

Cat. No. GTX78910

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-Fr
Reactivity	Human, Mouse, Rat, Monkey

Package
100 µg

Applications

Application Note

Recommended Starting Dilutions:

For IHC: Use at 1-3 µg/mL

Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GST-fusion protein aa260-378 of Presenilin 1.
Purification	Protein A/G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).