

uPA antibody [H77A10]

Cat. No. GTX79572

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	50mM Sodium Phosphate, 150mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length wild-type mouse protein (Glycosylated).
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX79572 WB Image

WB analysis of various sample lysates using GTX79572 uPA antibody [H77A10]

Lane 1 : Mouse uPA (0.5 µg) Non-reduced

Lane 2 : Mouse uPA (0.5 µg) + CPAI (2.5 µg) Non-reduced

Lane 3 : Prestained Standard

Dilution : 1:3000



For full product information, images and publications, please visit our [website](https://www.genetex.com).