

tPA antibody (FITC)

Cat. No. GTX79612

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Application	ELISA
Reactivity	Human, Pig

Package
1 mg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Not tested in other applications.	

Calculated MW 63 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length wild-type human protein (Glycosylated).
Purification	Protein G purified From serum
Conjugation	Fluorescein isothiocyanate (FITC)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).