

## SERPINE2 antibody (FITC)

## Cat. No. GTX79671

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Mouse

Package 500 μg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Properties		
Form	Liquid	
Buffer	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA	
Preservative	No preservatives	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.	
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)	
Immunogen	Recombinant full length wild-type mouse protein (Glycosylated).	
Purification	Protein G purified From serum	
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 1