

## Prorenin/Renin antibody (FITC)

**Cat. No. GTX79675**

<b>Host</b>	Sheep
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Application</b>	ELISA
<b>Reactivity</b>	Mouse, Rat

**Package**  
500 µg

### APPLICATION

#### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

<b>Calculated MW</b>	44 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	Recognises Renin/Prorenin

### PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA
<b>Preservative</b>	No preservative
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
<b>Concentration</b>	10 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant full length rat protein fused to a His tag at C-terminal end (Glycosylated).
<b>Purification</b>	IgG fraction From serum
<b>Conjugation</b>	Fluorescein isothiocyanate (FITC)
<b>Note</b>	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).