

Prorenin/Renin antibody (FITC)

Cat. No. GTX79675

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Mouse, Rat

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note Recognises Renin/Prorenin

Properties

Form	Liquid
Buffer	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length rat protein fused to a His tag at C-terminal end (Glycosylated).
Purification	IgG fraction From serum
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).